


# BE-DOUBLEFLOW



**SPECIFICATIONS**

**TECHNICAL FEATURES**

**FILLING AND HANDLING**


**COATING**

**CELL CULTURE**

**FLOW SET UP**

**CULTURE MODELS**

# BE-DOUBLEFLOW



This is our **most advanced** device. It consists of **two perfusable channels connected via a porous membrane**. Explore the crosstalk between different **2D and 3D cultures in a biomimetic environment** and **control** the efficiency of the interaction by selecting the **optimal pore size** for your application. This is the optimal device when a **hypoxia environment** is needed, for studying **the effect of circulating particles** (bacteria, immune system, circulating tumor cells) and for **endothelium/epithelium barrier** when flux plays a role in both sides of the coculture.

## TECHNICAL FEATURES



The design of the BE-Doubleflow consists of **2 channels connected through a porous membrane**. Flow is possible in both channels with the help of screw-like designed inlet and outlet wells by screwing them with connectors and tubes to a fluidic system. Evaporation reservoirs are next to the medium reservoirs to be filled with PBS or water during the incubation before closing the system with tubes.

## CONTENT

The box contains **10 BE-Doubleflow** chips that are **individually packaged**. Each chip undergoes a sterilization process before sending the chip to the final user. Chips should be stored at room temperature in **dry places** which are **not exposed to direct sunlight** at room temperature (**15-25°C**).

	Height	Width	Length	Total volume
<b>Each channel</b>	375 µm	1,5 mm	46 mm	31,2 µL
<b>Inlet/Outlet</b>	7 mm	UNF 1/4" - 28		130 µL
<b>Medium reservoir</b>	5 mm	3,6 mm	8,8 mm	185 µL
<b>Membrane pore size</b>	1 µm			

*\* The volumes presented on the table are theoretical values calculated for the standard products. Changes in the device features of custom chips may modify the exact channel volume.*



# BE-DOUBLEFLOW

## FILLING AND HANDLING

Before seeding, **prewarm the device** in the incubator overnight to avoid the appearance of air bubbles.

When filling the channels up, **place the tip of the pipette completely vertical** into the channel's inlet pinhole. Inject the liquid with a **continuous and constant flow** to avoid backward flow.



Use the **lateral holders** for easy and safe handling.



## COATING

The following protocol is applicable for channel **coating, washing, staining, cell extraction and medium refreshing**.

**1. Fill the channel volume up with the coating solution by pipetting through the pinhole and incubate at manufacturer's conditions.**





**Note:** Another way to proceed is to **pipette the solution into the inlet well** and let the liquid to flow by itself through the channel **until it reaches the outlet well**.

**Tip:** When seeding a monolayer underneath the membrane, flip the chip **downwards** and incubate for better coating deposition on the membrane. In this case the volume must be the one that fits in the channel.


**2. Remove the solution** from the reservoirs and inlet/outlet wells and **wash the channel** by adding the recommended dilution buffer into the inlet well and removing it at the outlet channel with the pipette. **Repeat** this step as many times as necessary.



**Tip:** You may **tilt the device** manually in order to make the dilution buffer flow **several times** from one reservoir to another before you remove it.



**3. Aspirate** the dilution buffer completely before seeding.



**Important:** Note that the channel should never dry up. Take incubation times into account and use the evaporation reservoirs as long as possible.

# BE-DOUBLEFLOW

## CELL SEEDING/CELL CULTURE

**1. Trypsinize and count cells** as usual. Cell concentration will vary with cell type and channel dimensions. It is recommended to seed a cell concentration to **obtain a confluent layer within 2-3 days**.

### Seeding underneath the membrane

**2. Fill the lower channel up** with the exact volume of the channel by pipetting through the pinhole. Find the **exact volume** of your chip in the table of the technical feature section.

**Important:** Adding the **correct volume** into the channel during seeding **helps attaching cells** to the lower side of the membrane.

**3.** Cover the inlets, **flip the chip downwards** and incubate until the cells attach (time may vary depending on the cell type).

**4.** After cell attachment, **add more medium** to the medium reservoirs. **Cover and incubate** until needed.

**3D culture:** In case of **seeding a hydrogel, flip the chip** upwards and downwards **every 10 seconds** for a few minutes (depending on the hydrogel polymerization time) for **homogeneous 3D distribution**.

**3D culture** is possible to perform both in upper and lower channel. BEOnChip recommends **seeding hydrogels into the lower channel** for better visualization of cells under inverted **microscopy**.

### Seeding on the membrane

**5. Fill the upper channel up** with the exact volume of the channel by pipetting through the pinhole. Find the **exact volume** of your chip in the table of the technical feature section.

**Important:** Adding the **correct volume** into the channel during seeding on the upper channel **helps saving cells** and makes possible to **concentrate them**. This will allow to create a **cell monolayer** along the channel easily/quickly.

**6.** Add **PBS or water** to the **evaporation reservoirs** and incubate until the cells are completely attached (time may vary depending on the cell type).

**7.** Add more **medium to the reservoirs, cover and incubate** until reaching confluence.

**8. Refresh medium** according to cellular requirements. We recommend to **repeat this process** every one or two days.

**Important:** Be aware **not to remove the medium inside of the channel** during medium refreshing. This will prevent you to **avoid the entrance of air bubbles**. Thus, when removing the solution, **place the tip of the pipette away from the pinhole** to avoid emptying the channel.

For optimal seeding inside the channel, we recommend a constant renewal of the medium once the cells have adhered to the surface of the device (from two to six hours after seeding).

## FLOW SET UP

Before setting the flow up:

- **Sterilize and prewarm** the tubes and the fluidic elements overnight at 37°C.
- Set the system in a **laminar flow cabinet**.
- **Cells must be well adhered** to the surface before mounting the perfusion system.
- The channels and inlet/outlet wells should **never be depleted of culture medium**.
- Both inlets and outlets are designed to be for **screw connectors (1/4" - 28)**.

For more information regarding [connecting a microfluidic flow system check our website](#).

### Connecting a perfusion system

**9. Remove the medium** from the reservoirs.


**10. Screw the outlet connector** and ensure that the tube is perfectly fixed. **Remove the displaced medium** from the reservoir.

**11. Prime the tube** that will **connect to the inlet** before assembling. This will **prevent air bubbles** from entering the device.

**12. Screw the inlet connector** and remove the displaced medium from the reservoirs.

**13.** Once the system is closed, **switch the flow on**.

**Check** that there are **no leaks** in the system by leaving the pump **running for a couple of minutes** before placing the device in the incubator. Medium reservoirs will be filled up if the connections are not completely closed.



**Important:** Note that the channel should never dry up. Take incubation times into account and use the evaporation reservoirs as long as possible.


# BE-DOUBLEFLOW

## CULTURE MODELS


Depending on the assay, either **monoculture or coculture can be performed**. Find below several of the **most frequently used culture models** on the BE-Doubleflow:

### Seeding underneath the membrane



**1) Seed a monolayer over and/or below** the membrane with different cell types.



**2) For more realistic assays it is possible to add a hydrogel for 3D cell culture or combine 2D and 3D.**



**3) Perfuse immune system cells, CTC's or bacteria** and explore their interaction with the culture.



**BEOnChip S.L.**  
Campus Río Ebro  
Mariano Esquillor s/n  
Edificio CEMINEM  
50018, Zaragoza, Spain

[info@beonchip.com](mailto:info@beonchip.com)  
+34 655 16 16 91

[www.beonchip.com](http://www.beonchip.com)