



Aria Software Development Kit

User Manual

Contents

1	Introduction to Aria SDK	7
2	Requirements	8
2.1	System requirements	8
2.2	Supported instrument	8
3	SDK general philosophy	9
3.1	Pipeline	9
3.2	Status management	9
4	Software layers	10
4.1	Aria SDK native shared libraries	10
4.2	Middleware	10
4.3	Installation	11
4.3.1	C++	11
4.3.2	C#	11
4.3.3	Python	11
5	Fluigent SDK Functions	13
5.1	Types definition	13
5.2	SDK Wrapper	15
5.2.1	GetErrorMessage	15
5.2.2	GetErrorSeverity	15
5.2.3	GetErrorTimestamp	15
5.2.4	ResetErrors	16
5.2.5	HasAsyncError	16
5.2.6	GetAsyncErrorCount	16
5.2.7	TryGetNextAsyncError	16

5.2.8	LoadPhysicalInstrument	17
5.2.9	LoadSimulatedInstrument	17
5.2.10	IsInstrumentSimulated	17
5.2.11	GetFlowUnitType	17
5.2.12	GetExternalSwitchType	18
5.2.13	GetAriaSerialNumber	18
5.2.14	GetFirmwareVersion	18
5.2.15	GetFlowEZFirmwareVersion	18
5.2.16	GetMinFlowRate	19
5.2.17	GetMaxFlowRate	19
5.2.18	GetMinPressure	19
5.2.19	GetMaxPressure	19
5.2.20	GetDateTimeFormat	20
5.2.21	SetPrefillAndPreloadFlowRate	20
5.2.22	GetPrefillAndPreloadFlowRatePreset	20
5.2.23	GetPrefillAndPreloadFlowRate	20
5.2.24	SetCalibrationValue	21
5.2.25	SetStep3CalibrationValue	21
5.2.26	SetCalibrationValues	21
5.2.27	GetCalibrationValues	21
5.2.28	GetStep3CalibrationValues	22
5.2.29	GetMaxStep3CalibrationValueCount	22
5.2.30	StartSequence	23
5.2.31	GenerateSequenceJSON	23
5.2.32	LoadSequenceFromJSON	23
5.2.33	LoadSequence	23
5.2.34	GetReservoirEstimatedRequiredVolume	24
5.2.35	IsReservoirEstimatedOverCapacity	24
5.2.36	EnablePrefill	24
5.2.37	IsPrefillEnabled	24
5.2.38	EnableZeroPessureMode	25
5.2.39	IsZeroPressureModeEnabled	25
5.2.40	SetSequenceStartASAP	25

5.2.41 IsSequenceStartingASAP	25
5.2.42 SetSequenceStartTime	26
5.2.43 GetSequenceStartTime	26
5.2.44 GetTotalDuration	26
5.2.45 GetSequenceStepCount	27
5.2.46 RemoveStep	27
5.2.47 InsertFlushStep	27
5.2.48 InsertSendSignalStep	28
5.2.49 InsertTimedInjectionStep	28
5.2.50 InsertVolumeInjectionStep	28
5.2.51 InsertWaitStep	29
5.2.52 InsertWaitUserStep	29
5.2.53 InsertWaitSignalStep	30
5.2.54 GetEstimatedStepStartTime	30
5.2.55 GetEstimatedStepDuration	31
5.2.56 GetStepType	31
5.2.57 SetParameter	31
5.2.58 SetParameter	31
5.2.59 SetParameter	32
5.2.60 SetParameter	32
5.2.61 SetParameter	32
5.2.62 GetIntParameter	33
5.2.63 GetBoolParameter	33
5.2.64 GetFloatParameter	33
5.2.65 GetStringParameter	34
5.2.66 GetSignalTypeParameter	34
5.2.67 SetFlowRateOrder	35
5.2.68 GetFlowRateOrder	35
5.2.69 GetMeasuredFlowRate	35
5.2.70 SetPressureOrder	35
5.2.71 GetPressureOrder	36
5.2.72 GetMeasuredPressure	36
5.2.73 SelectReservoir	36

5.2.74 GetSelectedReservoir	36
5.2.75 StopFlow	37
5.2.76 IsFlowStopped	37
5.2.77 GetExternalSwitchMaxReachablePort	37
5.2.78 SetExternalSwitchPort	37
5.2.79 GetCurrentExternalSwitchPort	38
5.2.80 SetEnabledPort	38
5.2.81 IsPortEnabled	38
5.2.82 GetWastePort	38
5.2.83 GetDefaultOutputPort	39
5.2.84 IsSequenceRunning	40
5.2.85 PauseSequence	40
5.2.86 IsSequencePaused	40
5.2.87 ResumeSequenceExecution	40
5.2.88 Cancel	41
5.2.89 GetPrefillStepNumber	41
5.2.90 GetPreloadStepNumber	41
5.2.91 GetCurrentStep	41
5.2.92 GetProgress	42
5.2.93 GetPrefillAndPreloadProgress	42
5.2.94 HasSequenceEnded	42
5.2.95 GetLastMeasuredCalibrationVolume	43
5.2.96 GetCalibrationState	43
5.2.97 StartCalibrationStep1	43
5.2.98 StartCalibrationStep2	43
5.2.99 StartCalibrationStep3_2Switch	44
5.2.100StartCalibrationStep3_MSwitch	44
5.2.101ValidateCalibration	44
5.2.102CancelCalibration	44
5.2.103StartCleaning1_Water	44
5.2.104StartCleaning2_Tergazyme	45
5.2.105StartCleaning3_Air	45
5.2.106StartCleaning4_IPA	45

5.2.107StartCleaning5_Air	45
5.2.108CancelCleaning	46
5.2.109SendTTLSignal	47
5.2.110StartAwaitingTTLSignal	47
5.2.111StartAwaitingTTLSignal	47
5.2.112CheckTTLSignal	47
5.2.113StopAwaitingTLL	48
5.2.114SetTTLPulseDuration	48
5.2.115GetTTLPulseDuration	48
5.2.116SendTCPMessage	48
5.2.117StartAwaitingTCPMessage	49
5.2.118CheckTCPMessage	49
5.2.119StopAwaitingTCPMessage	49
5.2.120SetTCPMode	49
5.2.121IsTCPServerMode	50
5.2.122SetTCPPort	50
5.2.123GetTCPPort	50

1 | Introduction to Aria SDK

Aria Software Development Kit (SDK) allows you to fully integrate Aria device in your application; it has been declined in several languages, namely C#.NET and, in beta version, C++ and Python.

The aim of this document is to introduce the SDK's exposed functions which can be used to interact with your Aria.

2 | Requirements

2.1 System requirements

The Aria SDK can only run on Windows systems for the moment. Any version more recent than Windows 10 (included) are supported.

2.2 Supported instrument

By using Aria SDK, you have direct access to Aria instrument as a whole, when using a sequence, or to certain of its individual components when used "remotely". The controllable components are:

- Flow EZ™ pressure controller
- FlowUnit M or L depending on the Aria model
- Internal M-Switch (reservoir selection)
- Internal 2-Switch (stop flow)
- External valve (M-Switch or 2-Switch depending on the Aria model)

3 | SDK general philosophy

As for the Aria UI software, the Aria SDK lies on the control of Aria via sequences (or protocols). Even if some remote control of the individual components is possible, the SDK really shines when it comes to schedule some long sequences of commands. Taking all possible parameters into account, Aria SDK is able to play a sequence of injection from multiple reservoirs and to multiple output channels while minimizing the consumption of chemical products (buffer, cell culture, etc.). The current manual aims to provide technical help to setup, control and monitor your Aria in a programmatic way. For any functioning details of the Aria instrument itself, scientific applications or hardware specificities, please refer to the Aria user manual.

3.1 Pipeline

An Aria experiment typically follows this suite of actions:

1. Instrument detection (real or simulated)
2. (opt.) Check configuration (SN, etc.)
3. Calibration (calculate the internal volumes for accurate prediction)
4. Sequence edition
5. (opt.) Save to file
6. Run sequence

3.2 Status management

When called, each function returns an error ID. If the command was properly executed a -1 value ID is returned, otherwise a new ID (incremented from the last error ID) is returned. All errors are saved in a stack during a session lifetime and none is cleared by default. Only a call to `ResetErrors` (see below) during a session would result in an empty stack.

The lone ID does not provide any details about the error returned (other than the presence of the error). To handle error details, three specific functions can be used:

- `GetErrorSeverity`: Returns the error severity as an `ErrorSeverity` enum.
- `GetErrorMessage`: Returns the error details as a string.
- `ResetErrors`: Cleans the error stack, next error ID will be 0.

4 | Software layers

The Aria SDK is based on a native library built for Windows and written in C#. This library handles low-level communication with the Aria instrument. Calling the native libraries directly is possible, but is recommended only for advanced users.

Additionally, more friendly packages and examples are provided for two major programming languages so far: C# and Python. They are collectively referred to as Middleware in this manual.

We strongly recommend using the Middleware if your programming language of choice is supported. It is open source, so you can modify it to suit your needs.

4.1 Aria SDK native shared libraries

The native shared libraries are provided in the **shared/** folder of the SDK archive, together with the C#.NET and Python examples, respectively in the **csharp/** and **python/** folders.

Any language that interfaces with C# should be able to access the library functions, as demonstrated in the Middleware source code.

The library functions are generally non-blocking and return immediately, with the exception of the functions that start a sequence or a procedure or that explicitly wait for a signal, such as `StartAwaitingTCPMessage` and `StartAwaitingTTLSignal`.

These values represent expected response time both when reading and when setting values on the instrument. Calling `GetXXX` functions more frequently than these delays will simply return the same value repeatedly until it is updated by the instrument. Calling `SetXXX` functions more frequently might cause the library to block while it waits for the instrument to process the commands.

For your information, The data refresh rate of Aria is **10ms**.

4.2 Middleware

The SDK middleware is a set of packages that make it easier to use the SDK with various programming languages. So far, they mainly act as examples and can be used as they are as a coding starting point for new developments.

The following programming languages are supported:

Language	Package
C++	[IN DEVELOPMENT]
C#	Program.cs example script aria-sdk-example.sln Visual Studio complete solution containing middleware and examples
Python	aria-sdk-example.py example script

The middleware matches the conventions of each programming language while keeping the interface as similar as possible across all supported languages.

The following sections contain installation and usage instructions for each language.

4.3 Installation

See for each language the installation specificities.

4.3.1 C++

[IN DEVELOPMENT]

4.3.2 C#

The C# middleware consists of a Visual Studio solution (**aria-sdk-example.sln**) containing:

- A .NET Framework 4.8 middleware **aria-sdk-example.csproj** project file
- An example script **Program.cs**
- A **README.md** file

Simply copy the DLL library to the example folder then build and run the project.

4.3.3 Python

The Python package groups:

- An example script **aria-sdk-example.py**
- A **README.md** file

You can easily take the first part of the example script (Initialization section) and add your code after it. Do not forget to change the .dll path if you move it away from the script location. Here is the line to edit:

```
path = os.path.join(os.getcwd(), 'aria-sdk.dll')
```

Language specifics:

- The Python support involves the usage of **pythonnet** package that can be found on PIP (<https://pypi.org/project/pythonnet/>).
- Functions that return values (such as the functions starting with GetXXX) returns both the desired value and the error as a tuple. You must store both value into a tuple or an exception will be raised.

```
stepProgress = Monitoring.GetProgress(currentStep)           # WRONG  
stepProgress, error = Monitoring.GetProgress(currentStep)    # OK
```

- For more details about how ref and out values are handled in **Python.NET**, please refer to <https://pythonnet.github.io/pythonnet/python.html#out-and-ref-parameters>

5 | Fluigent SDK Functions

5.1 Types definition

1. ErrorSeverity

Returned error severity when requested by GetErrorSeverity.

Value	Enum	Description
0	Info	No error
1	Warning	Report non-blocking warning
2	Error	Report failed action

2. FlowUnitType

Type of the internal FlowUnit. Only FlowUnits **M** and **L** are currently available for Aria.

Value	Enum	Description
0	UnknownFlowUnit	Cannot get the FlowUnit type
1	<i>XS</i>	<i>FlowUnit XS (NA)</i>
2	<i>S</i>	<i>FlowUnit S (NA)</i>
4	<i>M</i>	FlowUnit M
8	<i>L</i>	FlowUnit L
16	<i>XL</i>	<i>FlowUnit XL (NA)</i>
32	<i>MPLUS</i>	<i>FlowUnit M+ (NA)</i>
64	<i>LPLUS</i>	<i>FlowUnit L+ (NA)</i>

3. SignalType

Aria allows to send/receive TTL binary signals and TCP/IP messages. Those signals can be sent at the beginning and/or end of any sequence function. See Sequence edition section for more details.

Value	Enum	Description
0	TTL	TTL binary signal
1	TCP	TCP/IP message

4. FlowRatePreset

It is possible to tune the flowrate used for the prefill step of a sequence. The FlowRatePreset indicates the balance between precision and speed to be used in the SetPrefillAndPreloadFlowRate function.

Value	Enum	Description
0	Precision	FlowUnit M: 30 µl/min FlowUnit L: 50 µl/min
1	Balanced	FlowUnit M: 55 µl/min FlowUnit L: 250 µl/min
2	Fast	FlowUnit M: 80 µl/min FlowUnit L: 500 µl/min
3	Max	FlowUnit M: 80 µl/min FlowUnit L: 1000 µl/min

5. SwitchType

Type of the external Switch returned by GetExternalSwitchType.

Value	Enum	Description
0	UnknownSwitch	Cannot get the Switch type
1	TwoSwitch	2-Switch (3-port/2-way valve)
2	MSwitch	M-Switch (11-port/10-position valve)

6. StepType

Type of a sequence step.

Value	Enum	Description
0	Flush	Flushed the liquid remaining in the tubing to the waste
1	TimeInjection	Injection based on time and flowrate
2	VolumInjection	Injection based on volume and flowrate
3	Wait	Wait for a certain time
4	WaitForUser	Wait until user input
5	WaitForExternalSignal	Wait for an external signal of type SignalType before proceeding
6	SendExternalSignal	Send a signal of type SignalType

7. StepParameter

Type of a step parameter. Used in SetParameter and all Get*Parameter functions.

Value	Enum	Description
0	PRE_SIGNAL	bool
1	PRE_SIGNAL_TYPE	SignalType
2	POST_SIGNAL	bool
3	POST_SIGNAL_TYPE	SignalType
4	INPUT_RESERVOIR	int (1 -> 10)
5	OUTPUT_DESTINATION	int (1 -> 2 with 2-Switch, 1 -> 10 with M-Switch)
6	FLOWRATE	float
7	VOLUME	float
8	DURATION	int
9	SIGNAL_MESSAGE	string
10	AWAITED_SIGNAL_TYPE	SignalType
11	BACKTRACK	bool

8. CalibrationState

State of the Calibration phase.

Value	Enum	Description
0	NotRunning	No Calibration phase currently running
1	Flushing	Flushing in progress
2	SettingUp	Setting up in progress
3	Calibrating	Calibration in progress

5.2 SDK Wrapper

Errors

5.2.1 GetErrorMessage

```
string GetErrorMessage(int errorId);
```

Returns the message associated to the error with ID *errorId*

Parameters

errorId	int	Error ID
---------	-----	----------

Output

errorMsg	string	Error message
----------	--------	---------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.2 GetErrorSeverity

```
ErrorSeverity GetErrorSeverity(int errorId);
```

Returns the error severity (as ErrorSeverity) associated to the error with ID *errorId*

Parameters

errorSeverity	ErrorSeverity	Error severity
---------------	---------------	----------------

Output

errorMsg	string	Error message
----------	--------	---------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.3 GetErrorTimestamp

```
string GetErrorTimestamp(int errorId);
```

Returns the timestamp (with the following format: yyyy/MM/dd-HH:mm:ss) associated to the error with ID *errorId*

Parameters

errorSeverity	ErrorSeverity	Error severity
---------------	---------------	----------------

Output

errorMsg	string	Error message
----------	--------	---------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.4 ResetErrors

```
int ResetErrors(int *errorId);
```

Clears the error stack of all previous errors. Next error ID will be then 0.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.5 HasAsyncError

```
bool HasAsyncError(int *errorId);
```

Returns true if the error stack has one error or more.

Output

hasError	bool	Error stack has one error or more
----------	------	-----------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.6 GetAsyncErrorCount

```
int GetAsyncErrorCount(int *errorId);
```

Returns the number of errors in the error stack.

Output

nbError	int	Number of errors in the error stack
---------	-----	-------------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.7 TryGetNextAsyncError

```
int TryGetNextAsyncError(int *errorId);
```

Gets the ID of the first (oldest) error from the error stack.

Output

errorId	int	ID or the first (oldest) error of the error stack
---------	-----	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

Instrument

5.2.8 LoadPhysicalInstrument

```
int LoadPhysicalInstrument(FlowUnitType flowUnit, SwitchType externalSwitch);
```

Searches for a connected Aria instrument and loads it. Returns true if a connected Aria instrument was detected, false otherwise.

Output

success	bool	Has the instrument be successfully loaded or not.
---------	------	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.9 LoadSimulatedInstrument

```
int LoadSimulatedInstrument(FlowUnitType flowUnit, SwitchType externalSwitch);
```

Loads a Simulated Instrument with the given flowUnit and externalSwitch types.

Parameters

flowUnit	FlowUnitType	Type of the FlowUnit to be simulated
externalSwitch	SwitchType	Type of the external Switch to be simulated

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.10 IsInstrumentSimulated

```
int IsInstrumentSimulated(bool *isSimulated);
```

Reports if the current instrument is simulated.

Output

isSimulated	bool	Is the instrument simulated or not.
-------------	------	-------------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.11 GetFlowUnitType

```
int GetFlowUnitType(FlowUnitType *type);
```

Returns the current Instrument FlowUnit Type.

Output

type	FlowUnitType	Type of the current FlowUnit.
------	--------------	-------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.12 GetExternalSwitchType

```
int GetExternalSwitchType(SwitchType *type);
```

Returns the current instrument external Switch type.

Output

type	SwitchType	Type of the external Switch type.
------	------------	-----------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.13 GetAriaSerialNumber

```
int GetAriaSerialNumber(int *serialNumber);
```

Returns Aria instrument Serial Number.

Output

serialNumber	int	Aria instrument SN.
--------------	-----	---------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.14 GetFirmwareVersion

```
int GetFirmwareVersion(int *firmwareVersion);
```

Returns Aria instrument firmware version.

Output

firmwareVersion	int	Aria instrument firmware version.
-----------------	-----	-----------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.15 GetFlowEZFirmwareVersion

```
int GetFlowEZFirmwareVersion(bool *firmwareVersion);
```

Returns Aria FlowEZ firmware version.

Output

firmwareVersion	int	Aria FlowEZ firmware version.
-----------------	-----	-------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.16 GetMinFlowRate

```
int GetMinFlowRate(float *flowRate);
```

Returns the minimum flowrate order allowed by the current instrument.

Output

flowRate	float	Minimum flowrate possible with the current instrument.
----------	-------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.17 GetMaxFlowRate

```
int GetMaxFlowRate(float *flowRate);
```

Returns the maximum flowrate order allowed by the current instrument.

Output

flowRate	float	Maximum flowrate possible with the current instrument.
----------	-------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.18 GetMinPressure

```
int GetMinPressure(float *pressure);
```

Returns the minimum pressure order allowed by the current instrument.

Output

pressure	float	Minimum pressure possible with the current instrument.
----------	-------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.19 GetMaxPressure

```
int GetMaxPressure(float *pressure);
```

Returns the maximum pressure order allowed by the current instrument.

Output

pressure	float	Maximum pressure possible with the current instrument.
----------	-------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

Configuration

5.2.20 GetDateTimeFormat

```
int GetDateTimeFormat(string *timeFormat);
```

Returns the DateTime format used in Aria SDK functions.

Output

timeFormat	string	DateTime format used in SDK functions.
------------	--------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.21 SetPrefillAndPreloadFlowRate

```
int SetPrefillAndPreloadFlowRate(FlowratePreset flowratePreset);
```

Defines the prefill and preload flowrate from the given FlowratePreset. This impacts the precision vs speed balance for those two steps.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.22 GetPrefillAndPreloadFlowRatePreset

```
int GetPrefillAndPreloadFlowRatePreset(FlowratePreset *flowratePreset);
```

Returns the current prefill and preload flowrate preset as FlowratePreset.

Output

flowratePreset	FlowRatePreset	Current flowrate preset.
----------------	----------------	--------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.23 GetPrefillAndPreloadFlowRate

```
int GetPrefillAndPreloadFlowRate(float *flowrate);
```

Returns the current prefill and preload flowrate (in µl/min).

Output

flowrate	float	Current prefill and preload flowrate (in µl/min).
----------	-------	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.24 SetCalibrationValue

```
int SetCalibrationValue(int stepId, float volume);
```

Sets the internal volume **volume** (in μl) for step **stepId**.

Parameters

stepId	int	Step ID
volume	float	Internal volume for step <i>step</i>

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.25 SetStep3CalibrationValue

```
int SetStep3CalibrationValue(int step3PortId, float volume);
```

Sets the internal volume **volume** (in μl) for port **step3PortId** of step 3 (2-Switch: 1 -> 2, M-Switch: 1 -> 10).

Parameters

step3PortId	int	Step 3 port ID (2-Switch: 1 -> 2, M-Switch: 1 -> 10)
volume	float	Internal volume for step <i>step</i>

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.26 SetCalibrationValues

```
int SetCalibrationValues(float step1Volume, float step2Volume, float[]  
    step3Volumes);
```

Sets the internal volumes **step1Volume** (in μl), **step2Volume** (in μl) and **step3Volumes** (in μl) for all Calibration steps.

Parameters

stepId	int	Step ID
volume	float	Internal volume for step <i>step</i> (in μl)

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.27 GetCalibrationValues

```
int GetCalibrationValues(float *step1Volume, float *step2Volume);
```

Returns the internal volumes (in μl) for Calibration steps 1 and 2.

Returns

step1Volume	float	Internal volume of Calibration step 1 (in μl).
step2Volume	float	Internal volume of Calibration step 2 (in μl).
errorId	int	Error ID (-1 if none)

5.2.28 GetStep3CalibrationValues

```
float[] GetStep3CalibrationValues(int *errorId);
```

Returns the internal volumes (in μl) for Calibration steps 1 and 2.

Output

step3Volumes	float[]	Internal volumes of Calibration step 3 (in μl).
--------------	---------	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.29 GetMaxStep3CalibrationValueCount

```
int GetMaxStep3CalibrationValueCount(int *nbMaxValues);
```

Returns the maximum number of calibration values in the Calibration step 3 table.

Output

nbMaxValues	int	Maximum number of calibration values in Calibration step 3 table.
-------------	-----	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

Sequence configuration

5.2.30 StartSequence

```
int StartSequence(int *errorId);
```

Start the current sequence.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.31 GenerateSequenceJSON

```
string GenerateSequenceJSON(int *errorId);
```

Returns the sequence saved as a JSON string. **NOT AVAILABLE IN PYTHON**

Output

sequenceAsJSON	string	JSON string representation of the current sequence
----------------	--------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.32 LoadSequenceFromJSON

```
int LoadSequenceFromJSON(string jsonString);
```

Load a sequence from a JSON String. **NOT AVAILABLE IN PYTHON**

Parameter

jsonString	string	JSON string representation of the current sequence
------------	--------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.33 LoadSequence

```
int LoadSequence(string filePath);
```

Load sequence from the file at **filePath**. **NOT AVAILABLE IN PYTHON**

Parameter

filePath	string	Path of a JSON file containing a sequence information
----------	--------	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.34 GetReservoirEstimatedRequiredVolume

```
float GetReservoirEstimatedRequiredVolume(int reservoirNumber);
```

Returns the estimated required volume with which to fill the given reservoir before starting the sequence.

Parameter

reservoirNumber	int	Reservoir ID (1 -> 10) to be considered.
-----------------	-----	--

Output

requiredVolume	float	Estimated required volume for the reservoir considered.
----------------	-------	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.35 IsReservoirEstimatedOverCapacity

```
bool IsReservoirEstimatedOverCapacity(int reservoirNumber);
```

Returns true if the required volume for the given reservoir is above that reservoir capacity, false otherwise. Reservoirs over capacity will require to refill them during the sequence execution.

Parameter

reservoirNumber	int	Reservoir ID (1 -> 10) to be considered.
-----------------	-----	--

Output

overCapacity	bool	Is the reservoir over capacity?
--------------	------	---------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.36 EnablePrefill

```
int EnablePrefill(bool enabled);
```

Enables or disables the prefill phase at the beginning of the sequence. The prefill fills the tubing between the reservoir and the internal M-Switch. Prefill can be safely disabled only if the tubing is already filled with the correct content.

Parameter

enabled	bool	Enables or not the prefill.
---------	------	-----------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.37 IsPrefillEnabled

```
bool IsPrefillEnabled(int *errorId);
```


Returns true if Prefill is enabled, false otherwise.

Output

enabled	bool	Is the prefill enabled (true) or not (false)
---------	------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.38 EnableZeroPessureMode

```
int EnableZeroPessureMode(bool enabled);
```

Enables or disables Zero Pressure mode. Zero Pressure mode forces the pressure to reset to 0 every time a Switch is moved to avoid flow rate spikes, irregularities, etc.

Parameter

enabled	bool	Enables or not the Zero Pressure mode.
---------	------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.39 IsZeroPressureModeEnabled

```
int IsZeroPressureModeEnabled(int *errorId);
```

Returns true if Zero Pressure Mode is enabled, false otherwise.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.40 SetSequenceStartASAP

```
int SetSequenceStartASAP(bool startASAP);
```

Defines if the sequence must start as soon as possible, or with a delay.

Parameter

startASAP	bool	Start the sequence as soon as possible or not.
-----------	------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.41 IsSequenceStartingASAP

```
bool IsSequenceStartingASAP(int *errorId);
```

Returns true if *startASAP* is enabled, false otherwise.

Output

startingASAP	bool	Will the sequence start ASAP?
--------------	------	-------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.42 SetSequenceStartTime

```
int SetSequenceStartTime(string dateTime);
```

Defines the time at which the sequence will be executed if *startASAP* is disabled. Time must be in the following format: yyyy/MM/dd-HH:mm:ss

Parameter

dateTime	string	Time at which the sequence will be started.
----------	--------	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.43 GetSequenceStartTime

```
string GetSequenceStartTime(int *errorId);
```

Returns the estimated start time of the first step of the sequence.

Output

dateTime	string	Time at which the sequence will be started.
----------	--------	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.44 GetTotalDuration

```
int GetTotalDuration(int *errorId);
```

Returns the estimated total duration of the sequence.

Output

duration	int	Sequence estimated total duration (in seconds).
----------	-----	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

Sequence edition

5.2.45 GetSequenceStepCount

```
int GetSequenceStepCount(int *errorId);
```

Returns the total number of steps in the current sequence.

Output

nbSteps	int	Number of current sequence steps
---------	-----	----------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.46 RemoveStep

```
int RemoveStep(int index);
```

Removes the step at **index** from the sequence.

Parameter

index	int	Index of the step to be removed
-------	-----	---------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.47 InsertFlushStep

```
int InsertFlushStep(int index, int inputReservoir, float flowRate, bool preSignal, SignalType preSignalType, bool postSignal, SignalType postSignalType);
```

Inserts a step of flushing in the current sequence at **index**. Reservoir **inputReservoir** will be flushed at **flowRate** µl/min.

Parameters

index	int	Position in the sequence where the step will be inserted. 0 to insert it at the beginning, -1 to insert it at the end.
inputReservoir	int	Input reservoir to be used (1 -> 10).
flowRate	float	Flowrate order to reach for this step (in µl/min).
preSignal	bool	Send a pre-signal.
preSignalType	SignalType	Pre-signal type.
postSignal	bool	Send a post-signal.
postSignalType	SignalType	Post-signal type.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.48 InsertSendSignalStep

```
int InsertSendSignalStep(int index, string message, bool preSignal, SignalType preSignalType, bool postSignal, SignalType postSignalType);
```

Inserts a step of TCP/IP signal sending (with message **message**) in the current sequence at **index**.

Parameters

index	int	Position in the sequence where the step will be inserted. 0 to insert it at the beginning, -1 to insert it at the end.
message	string	TCP/IP message to be sent.
preSignal	bool	Send a pre-signal.
preSignalType	SignalType	Pre-signal type.
postSignal	bool	Send a post-signal.
postSignalType	SignalType	Post-signal type.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.49 InsertTimedInjectionStep

```
int InsertTimedInjectionStep(int index, int inputReservoir, int destination, float flowRate, int duration_s, bool preSignal, SignalType preSignalType, bool postSignal, SignalType postSignalType);
```

Inserts a step of timed injection in the current sequence at **index**. Input reservoir **inputReservoir** will be injected into output port **destination** at **flowRate** µl/min for **duration_s** seconds.

Parameters

index	int	Position in the sequence where the step will be inserted. 0 to insert it at the beginning, -1 to insert it at the end.
inputReservoir	int	Input reservoir to be used (1 -> 10).
destination	int	External Switch port to be used (1 -> 2 for 2-Switch, 1 -> 10 for M-Switch).
flowRate	float	Flowrate order to reach for this step (in µl/min).
duration_s	int	Time of the injection (in seconds).
preSignal	bool	Send a pre-signal.
preSignalType	SignalType	Pre-signal type.
postSignal	bool	Send a post-signal.
postSignalType	SignalType	Post-signal type.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.50 InsertVolumeInjectionStep

```
int InsertVolumeInjectionStep(int index, int inputReservoir, int destination, float flowRate, int volume, bool preSignal, SignalType preSignalType, bool postSignal, SignalType postSignalType);
```

Inserts a step of volume injection in the current sequence at **index**. A volume of **volume** µl will be injected from input reservoir **inputReservoir** to output port **destination** at **flowRate** µl/min.

Parameters

index	int	Position in the sequence where the step will be inserted. 0 to insert it at the beginning, -1 to insert it at the end.
inputReservoir	int	Input reservoir to be used (1 -> 10).
destination	int	External Switch port to be used (1 -> 2 for 2-Switch, 1 -> 10 for M-Switch).
flowRate	float	Flowrate order to reach for this step (in µl/min).
volume	float	Volume to be injected (in µl).
preSignal	bool	Send a pre-signal.
preSignalType	SignalType	Pre-signal type.
postSignal	bool	Send a post-signal.
postSignalType	SignalType	Post-signal type.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.51 InsertWaitStep

```
int InsertFlushStep(int index, int duration_s, bool preSignal,
    SignalType preSignalType, bool postSignal, SignalType
    postSignalType);
```

Inserts a step of waiting in the current sequence at **index**. Step will wait for **duration_s** seconds before proceeding to the next step.

Parameters

index	int	Position in the sequence where the step will be inserted. 0 to insert it at the beginning, -1 to insert it at the end.
duration_s	int	Waiting time (in seconds).
preSignal	bool	Send a pre-signal.
preSignalType	SignalType	Pre-signal type.
postSignal	bool	Send a post-signal.
postSignalType	SignalType	Post-signal type.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.52 InsertWaitUserStep

```
int InsertWaitUserStep(int index, int timeout_s, bool preSignal, SignalType
    preSignalType, bool postSignal, SignalType postSignalType);
```

Inserts a step of waiting for user in the current sequence at **index**. Step will wait for the execution of ResumeSequenceExecution before proceeding to the next step.

Parameters

index	int	Position in the sequence where the step will be inserted. 0 to insert it at the beginning, -1 to insert it at the end.
timeout_s	int	Waiting timeout if no signal has been received (in seconds).
preSignal	bool	Send a pre-signal.
preSignalType	SignalType	Pre-signal type.
postSignal	bool	Send a post-signal.
postSignalType	SignalType	Post-signal type.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.53 InsertWaitSignalStep

```
int InsertWaitSignalStep(int index, int timeout_s, SignalType signalType, bool
    enableBacktrack, bool preSignal, SignalType preSignalType, bool postSignal,
    SignalType postSignalType);
```

Inserts a step of waiting for signal in the current sequence at **index**. Step will wait for a signal of type **SignalType** before proceeding to the next step. When **enableBacktrack** is true

Parameters

index	int	Position in the sequence where the step will be inserted. 0 to insert it at the beginning, -1 to insert it at the end.
timeout_s	int	Waiting timeout if no signal has been received (in seconds).
signalType	SignalType	Type of the signal to be waiting for.
preSignal	bool	Send a pre-signal.
preSignalType	SignalType	Pre-signal type.
postSignal	bool	Send a post-signal.
postSignalType	SignalType	Post-signal type.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.54 GetEstimatedStepStartTime

```
int GetEstimatedStepStartTime(int index);
```

Returns the estimated delay (in seconds) before the step at **index** is executed.

Parameter

index	int	Index of the selected step.
-------	-----	-----------------------------

Output

delay	int	Number of seconds before the selected step will be executed.
-------	-----	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.55 GetEstimatedStepDuration

```
int GetEstimatedStepDuration(int index);
```

Returns the estimated duration (in seconds) of the Step at **index**.

Parameter

index	int	Index of the selected step.
-------	-----	-----------------------------

Output

duration	int	Estimated duration of the selected step (in seconds).
----------	-----	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.56 GetStepType

```
StepType GetStepType(int index);
```

Returns the Step type (as StepType) of the step at **index**.

Parameter

index	int	Index of the selected step.
-------	-----	-----------------------------

Output

stepType	StepType	Step type of the selected step.
----------	----------	---------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.57 SetParameter

```
int SetParameter(int index, StepParameter parameterType, int value);
```

Assigns value **value** of the parameter **parameterType** (StepParameter) of the step at **index**. Check that the value type corresponds to the type of the step parameter (see StepParameter).

Parameters

index	int	Index of the selected step.
parameterType	StepParameter	Type of the step parameter to edit.
value	int	New value of the step parameter.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.58 SetParameter

```
int SetParameter(int index, StepParameter parameterType, bool value);
```

Assigns value **value** of the parameter **parameterType** (StepParameter) of the step at **index**. Check that the value type corresponds to the type of the step parameter (see StepParameter).

Parameters

index	int	Index of the selected step.
parameterType	StepParameter	Type of the step parameter to edit.
value	bool	New value of the step parameter.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.59 SetParameter

```
int SetParameter(int index, StepParameter parameterType, float value);
```

Assigns value **value** of the parameter **parameterType** (StepParameter) of the step at **index**. Check that the value type corresponds to the type of the step parameter (see StepParameter).

Parameters

index	int	Index of the selected step.
parameterType	StepParameter	Type of the step parameter to edit.
value	float	New value of the step parameter.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.60 SetParameter

```
int SetParameter(int index, StepParameter parameterType, string value);
```

Assigns value **value** of the parameter **parameterType** (StepParameter) of the step at **index**. Check that the value type corresponds to the type of the step parameter (see StepParameter).

Parameters

index	int	Index of the selected step.
parameterType	StepParameter	Type of the step parameter to edit.
value	string	New value of the step parameter.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.61 SetParameter

```
int SetParameter(int index, StepParameter parameterType, SignalType value);
```

Assigns value **value** of the parameter **parameterType** (StepParameter) of the step at **index**. One must check that the value type corresponds to the type of the step parameter (see StepParameter).

Parameters

index	int	Index of the selected step.
parameterType	StepParameter	Type of the step parameter to edit.
value	SignalType	New value of the step parameter.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.62 GetIntParameter

```
int GetIntParameter(int index, StepParameter parameterType);
```

Returns the parameter value (int) according to the **parameterType** for the step at **index**.

Parameters

index	int	Index of the selected step.
parameterType	StepParameter	Type of the step parameter.

Output

parameterValue	int	Value of the selected parameter.
----------------	-----	----------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.63 GetBoolParameter

```
bool GetBoolParameter(int index, StepParameter parameterType);
```

Returns the parameter value (bool) according to the **parameterType** for the step at **index**.

Parameters

index	int	Index of the selected step.
parameterType	StepParameter	Type of the step parameter.

Output

parameterValue	bool	Value of the selected parameter.
----------------	------	----------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.64 GetFloatParameter

```
float GetFloatParameter(int index, StepParameter parameterType);
```

Returns the parameter value (float) according to the **parameterType** for the step at **index**.

Parameters

index	int	Index of the selected step.
parameterType	StepParameter	Type of the step parameter.

Output

parameterValue	float	Value of the selected parameter.
----------------	-------	----------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.65 GetStringParameter

```
string GetStringParameter(int index, StepParameter parameterType);
```

Returns the parameter value (string) according to the **parameterType** for the step at **index**.

Parameters

index	int	Index of the selected step.
parameterType	StepParameter	Type of the step parameter.

Output

parameterValue	string	Value of the selected parameter.
----------------	--------	----------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.66 GetSignalTypeParameter

```
SignalType GetSignalTypeParameter(int index, StepParameter parameterType);
```

Returns the parameter value (SignalType) according to the **parameterType** for the step at **index**.

Parameters

index	int	Index of the selected step.
parameterType	StepParameter	Type of the step parameter.

Output

parameterValue	SignalType	Value of the selected parameter.
----------------	------------	----------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

Direct control

5.2.67 SetFlowRateOrder

```
int SetFlowRateOrder(float flowrate);
```

Set the flowrate order of the Aria instrument to **flowrate** (in $\mu\text{l}/\text{min}$).

Parameter

flowrate	float	Flowrate order to be reached (in $\mu\text{l}/\text{min}$)
----------	-------	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.68 GetFlowRateOrder

```
int GetFlowRateOrder(int *errorId);
```

Returns the last flowrate order sent to the Aria instrument (in $\mu\text{l}/\text{min}$).

Output

flowrate	float	Last flowrate order sent (in $\mu\text{l}/\text{min}$)
----------	-------	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.69 GetMeasuredFlowRate

```
int GetMeasuredFlowRate(int *errorId);
```

Returns the current flowrate value measured by the Aria instrument (in $\mu\text{l}/\text{min}$).

Output

flowrate	float	Current flowrate (in $\mu\text{l}/\text{min}$)
----------	-------	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.70 SetPressureOrder

```
int SetPressureOrder(float pressure);
```

Set the pressure order of the Aria instrument to **pressure** (in mBar).

Parameter

pressure	float	Pressure order to be reached (in mBar)
----------	-------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.71 GetPressureOrder

```
int GetPressureOrder(int *errorId);
```

Returns the last pressure order sent to the Aria instrument (in mBar).

Output

pressure	float	Last pressure order sent (in mBar)
----------	-------	------------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.72 GetMeasuredPressure

```
int GetMeasuredPressure(int *errorId);
```

Returns the current pressure value measured by the Aria instrument (in mBar).

Output

pressure	float	Current flowrate (in mBar)
----------	-------	----------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.73 SelectReservoir

```
int SelectReservoir(int reservoirId);
```

(Remote Control) Switch the internal M-Switch to connect to the given **reservoirId** (1 -> 10).

Parameter

reservoirId	int	Selected reservoir ID (1 -> 10)
-------------	-----	---------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.74 GetSelectedReservoir

```
int GetSelectedReservoir(int *errorId);
```

Returns the current selected reservoir.

Output

reservoirId	int	ID of the current selected reservoir
-------------	-----	--------------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.75 StopFlow

```
int StopFlow(bool stop);
```

(Remote Control) Switch the Internal 2-Switch to stop (true) or allow (false) flow.

Parameter

stop	bool	Stop (true) or allow (false) flow.
------	------	------------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.76 IsFlowStopped

```
bool IsFlowStopped(int *errorId);
```

Returns true if the flow is stopped by the Internal 2-Switch, false if it is open.

Output

stopped	bool	Is flow stopped (true) or not (false)
---------	------	---------------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.77 GetExternalSwitchMaxReachablePort

```
int GetExternalSwitchMaxReachablePort(int *errorId);
```

Returns the maximum Port number of the external Switch (2-Switch: 2 | M-Switch: 10)

Output

maxNumberPorts	int	Maximum number of ports for the current external Switch
----------------	-----	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.78 SetExternalSwitchPort

```
int SetExternalSwitchPort(int port);
```

(Remote Control) Switch the external Switch to the given chip **port**.

Parameter

port	int	Selected chip (external Switch) port.
------	-----	---------------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.79 GetCurrentExternalSwitchPort

```
int GetCurrentExternalSwitchPort(int *errorId);
```

Returns the current chip port of the external Switch.

Output

portId	int	ID of the current selected chip port (external Switch port)
--------	-----	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.80 SetEnabledPort

```
int SetEnabledPort(int port, bool enabled);
```

Enables or disables the given external Switch **port**. Enabled ports can be used in sequences as well as Calibration and Cleaning procedures.

Parameters

port	int	Selected chip (external Switch) port.
enabled	bool	Enables (true) or disables (false) selected port.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.81 IsPortEnabled

```
bool IsPortEnabled(int portId);
```

Returns true if the external Switch **port** is enabled, false otherwise.

Output

enabled	bool	Is external Switch port enabled (true) or not (false)
---------	------	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.82 GetWastePort

```
int GetWastePort(int *errorId);
```

Returns the external Switchp port number of the Waste port. (2-Switch: 2 | M-Switch: 10).

Output

portId	int	ID of the external Switch port used for the Waste.
--------	-----	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.83 GetDefaultOutputPort

```
int GetDefaultOutputPort(int *errorId);
```

Returns the external Switchp port number of the default Output port. (2-Switch: 1 | M-Switch: 1).

Output

portId	int	ID of the external Switch port used for the default Output.
--------	-----	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

Sequence monitoring

5.2.84 IsSequenceRunning

```
bool IsSequenceRunning(int *errorId);
```

Returns true if the current sequence is in progress, false otherwise. A paused sequence is still considered as running.

Output

running	bool	Is the sequence running (true) or not (false)
---------	------	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.85 PauseSequence

```
int PauseSequence(int *errorId);
```

Pauses the current sequence execution.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.86 IsSequencePaused

```
bool IsSequencePaused(int *errorId);
```

Returns true if the sequence is paused, false otherwise.

Output

paused	bool	Is the sequence paused (true) or not (false)
--------	------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.87 ResumeSequenceExecution

```
int ResumeSequenceExecution(int *errorId);
```

Resumes the execution of a paused sequence.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.88 Cancel

```
int Cancel(int *errorId);
```

Cancels the current procedure or sequence.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.89 GetPrefillStepNumber

```
int GetPrefillStepNumber(int *errorId);
```

Returns the step index of the Prefill phase of the sequence (base 1).

Output

stepId	int	Step ID of the Prefill phase for the current sequence.
--------	-----	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.90 GetPreloadStepNumber

```
int GetPreloadStepNumber(int *errorId);
```

Returns the step index of the Preload phase of the sequence (base 1).

Output

stepId	int	Step ID of the Preload phase for the current sequence.
--------	-----	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.91 GetCurrentStep

```
int GetCurrentStep(int *errorId);
```

Returns the step index of the current step of the sequence (base 1).

Output

stepId	int	Step ID of the current step of the running sequence
--------	-----	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.92 GetProgress

```
float GetProgress(int stepId);
```

Returns the progress level (%) of the step at index **stepId**. (base 1).

Parameter

stepId	int	Step to be considered.
--------	-----	------------------------

Output

progress	float	Progress level of the selected step
----------	-------	-------------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.93 GetPrefillAndPreloadProgress

```
float GetPrefillAndPreloadProgress(int *errorId);
```

Returns the cumulated progress for the Prefill and Preload phases.

Output

progress	float	Cumulated progress for the Prefill and Preload phases
----------	-------	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.94 HasSequenceEnded

```
bool HasSequenceEnded(int *errorId);
```

Returns true if the sequence execution ended.

Output

ended	bool	Has the current sequence ended (true) or not (false)
-------	------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

Procedures

5.2.95 GetLastMeasuredCalibrationVolume

```
float GetLastMeasuredCalibrationVolume(int *errorId);
```

Returns the last internal volume calculated during the Calibration phase.

Output

volume	float	Last internal volume calculated during the Calibration phase
--------	-------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.96 GetCalibrationState

```
CalibrationState GetCalibrationState(int *errorId);
```

Returns the current state of the Calibration phase (as CalibrationState).

Output

state	CalibrationState	Current state of the Calibration phase
-------	------------------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.97 StartCalibrationStep1

```
int StartCalibrationStep1(int *errorId);
```

Starts the 1st step of the Calibration phase.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.98 StartCalibrationStep2

```
int StartCalibrationStep2(int *errorId);
```

Starts the 2nd step of the Calibration phase.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.99 StartCalibrationStep3_2Switch

```
int StartCalibrationStep3_2Switch(int *errorId);
```

Starts the 3rd step of the Calibration phase for Aria instrument with external 2-Switch.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.100 StartCalibrationStep3_MSwitch

```
int StartCalibrationStep3_MSwitch(int *errorId);
```

Starts the 3rd step of the Calibration phase for Aria instrument with external M-Switch.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.101 ValidateCalibration

```
int ValidateCalibration(int *errorId);
```

Validates the current Calibration phase (results in saving the calibration volumes calculated and stopping the phase).

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.102 CancelCalibration

```
int CancelCalibration(int *errorId);
```

Cancels the current Calibration phase (results in discarding the calibration volumes calculated and stopping the phase).

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.103 StartCleaning1_Water

```
int StartCleaning1_Water(int[] reservoirNumbers, int bufferReservoirNumber);
```

Starts the step 1 of the Cleaning procedure : the cleaned **reservoirNumbers** should be filled with WATER, as well as the **bufferReservoirNumber** (9 - 10).

Parameters

reservoirNumbers	int[]	IDs of the reservoirs to be cleaned
bufferReservoirNumber	int	Error ID (-1 if none)

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.104 StartCleaning2_Tergazyme

```
int StartCleaning2_Tergazyme(int[] reservoirNumbers, int bufferReservoirNumber);
```

Starts the step 2 of the Cleaning procedure : the cleaned **reservoirNumbers** should be filled with TERGAZYME, as well as the **bufferReservoirNumber** (9 - 10).

Parameters

reservoirNumbers	int[]	IDs of the reservoirs to be cleaned
bufferReservoirNumber	int	Error ID (-1 if none)

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.105 StartCleaning3_Air

```
int StartCleaning3_Air(int[] reservoirNumbers, int bufferReservoirNumber);
```

Starts the step 3 of the Cleaning procedure : the cleaned **reservoirNumbers** should be EMPTY, as well as the **bufferReservoirNumber** (9 - 10).

Parameters

reservoirNumbers	int[]	IDs of the reservoirs to be cleaned
bufferReservoirNumber	int	Error ID (-1 if none)

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.106 StartCleaning4_IPA

```
int StartCleaning4_IPA(int[] reservoirNumbers, int bufferReservoirNumber);
```

Starts the step 4 of the Cleaning procedure : the cleaned **reservoirNumbers** should be filled with IPA, as well as the **bufferReservoirNumber** (9 - 10).

Parameters

reservoirNumbers	int[]	IDs of the reservoirs to be cleaned
bufferReservoirNumber	int	Error ID (-1 if none)

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.107 StartCleaning5_Air

```
int StartCleaning5_Air(int[] reservoirNumbers, int bufferReservoirNumber);
```

Starts the step 5 of the Cleaning procedure : the cleaned **reservoirNumbers** should be EMPTY, as well as the **bufferReservoirNumber** (9 - 10).

Parameters

reservoirNumbers	int[]	IDs of the reservoirs to be cleaned
bufferReservoirNumber	int	Error ID (-1 if none)

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.108 CancelCleaning

```
int CancelCleaning(int *errorId);
```

Cancels the current Cleaning phase.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

External communication

5.2.109 SendTTLSignal

```
int SendTTLSignal(int *errorId);
```

Sends a TTL signal through the Aria instrument.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.110 StartAwaitingTTLSignal

```
int StartAwaitingTTLSignal(int period);
```

Starts waiting for a TTL signal. TTL check is done every **period** milliseconds (adapt with respect to incoming TTL pulse duration). Must be stopped by StopAwaitingTTL

Parameters

period	int	Check frequency (in ms)
--------	-----	-------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.111 StartAwaitingTTLSignal

```
int StartAwaitingTTLSignal(int period, int timeout);
```

Starts waiting for a TTL signal. TTL check is done every **period** milliseconds (adapt with respect to incoming TTL pulse duration). Stops awaiting after **timeout** milliseconds.

Parameters

period	int	Check frequency (in ms)
timeout	int	Timeout duration (in ms)

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.112 CheckTTLSignal

```
bool CheckTTLSignal(int *errorId);
```

Checks if a TTL signal was received. If it was, stops the current TTL waiting process.

Output

signalReceived	bool	Has a signal been received (true) or not (false)
----------------	------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.113 StopAwaitingTLL

```
int StopAwaitingTLL(int *errorId);
```

Stops the current TTL waiting process.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.114 SetTTLPulseDuration

```
int SetTTLPulseDuration(int duration);
```

Defines the duration of the TTL pulse sent by the Aria instrument (TTL pulse = **duration** * 100ms).

Parameter

duration	int	Number of 100ms periods for a single TTL pulse signal
----------	-----	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.115 GetTTLPulseDuration

```
int GetTTLPulseDuration(int *errorId);
```

Gets the current duration of the TTL pulse sent by the Aria instrument. (TTL pulse = **duration** * 100ms).

Output

duration	int	Number of 100ms periods for a single TTL pulse signal
----------	-----	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.116 SendTCPMessage

```
int SendTCPMessage(string message);
```

Sends a TCP text message.

Parameter

message	string	Message to be sent
---------	--------	--------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.117 StartAwaitingTCPMessage

```
int StartAwaitingTCPMessage(string awaitedMessage);
```

Start waiting for a TCP message with the content **awaitedMessage**.

Parameter

awaitedMessage	string	Content of the message to be waited for
----------------	--------	---

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.118 CheckTCPMessage

```
bool CheckTCPMessage(string awaitedMessage);
```

Checks if a TCP message with the content **awaitedMessage** was received. If it was, stops the process waiting for this message.

Parameter

awaitedMessage	string	Content of the message to be waited for
----------------	--------	---

Output

signalReceived	bool	Has a signal been received (true) or not (false)
----------------	------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.119 StopAwaitingTCPMessage

```
int StopAwaitingTCPMessage(string awaitedMessage);
```

Stops waiting for a TCP message with the content **awaitedMessage**.

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.120 SetTCPMode

```
int SetTCPMode(bool enableServer);
```

Stops the current TTL waiting process.

Parameter

enableServer	bool	Enables server (true) or client (false) mode
--------------	------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.121 IsTCPServerMode

```
bool IsTCPServerMode(int *errorId);
```

Checks if Aria is in TCP Server mode.

Output

serverMode	bool	Is Aria in TCP server mode (true) or not (false)
------------	------	--

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.122 SetTCPPort

```
int SetTCPPort(int port);
```

Define the TCP port used for the TCP client and server.

Parameter

port	int	TCP client and server port number
------	-----	-----------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------

5.2.123 GetTCPPort

```
int GetTCPPort(int *errorId);
```

Stops the current TTL waiting process.

Output

port	int	TCP client and server port number
------	-----	-----------------------------------

Returns

errorId	int	Error ID (-1 if none)
---------	-----	-----------------------



FLUIGENT

O'kabé bureaux

57-77, avenue de Fontainebleau

94270 Le Kremlin-Bicêtre

FRANCE

Phone: +331 77 01 82 68

Fax: +331 77 01 82 70

www.fluigent.com

Technical support:

support@fluigent.com

Phone : +331 77 01 82 65

General information:

contact@fluigent.com