



L-X Specifications

P/N: IVLSW1

Version : 1.0

Performance	
Internal volume	660 nL
Dead volume	None
Switching time	100 ms
Mechanical	
Dimensions	7x9x15 cm
Weight	475 g
Maximum Pressure	7 bar (100 psi)
Fittings	Flangeless (1/16" OD)
Gas compatibility	Dry, oil-free gas, air, any non corrosive or non explosive gas
Liquid compatibility	Aqueous solvent, oil, organic solvent, biological sample
Wetted materials	PEEK
Electrical	
Port communication	RJ45
Power supply	RJ45
Software control	OxyGen or SDK