



RX Specifications

P/N: ISRX21

Version : 1.0

Performance	
Minimum pressure	500 mbar
Maximum pressure	2500 mbar
Free flow (flow rate at 0 mbar output)	5.6 L/min
Flow rate at 2000 mbar	0.8 L/min
Mechanical	
Dimensions	106 x 96.5 x 61.5 mm
Weight	823 g
Air quality*	Dried and filtered
Electrical	
Power voltage	24 V
Maximum current	800 mA
Maximum power consumption	19 W

*the RX has a built in drying component, it will reduce the humidity of the air supply to prevent condensation. This component is designed to perform in room temperature, for more detailed information on air supply conditions of your setup please contact us.