


Cable outlet







Version industrielle

Hole for M3 screws

Fittings XP235

R18,00

$\phi 3,20$



CREE / CREATED	DATE	ECHELLE	TITRE / TITLE	
LLE	10/06/2022	1:2	FGT_DRWEST_FU_INDUS	
MATERIAU / MATERIAL		SYSTEME	REFERENCE	FINITION / FINISH COULEUR / COLOR
			_FLG_01688	

Okabe bureaux, 47 avenue de Fontainebleau
94270 Le Kremlin-Bicêtre | France



TOLERANCES / TOLERANCES			
FORMAT	LANGUE	PROPRIETE	PAGE
A4		All proprietary rights reserved by Fluigent. Unauthorized reproduction will be legally pursued.	1 / 1