







4 x M5x0.8 - 6H ↧ 5,00
⌵ ⌀ 5,05 X 90°, TOP SIDE
PNEUMATIC IN/OUTs





Indice	Date	Modification	Approuvé
A	10/02/2022	DEFINITION	LLE

Table de Révision


CREE / *CREATED*	DATE	ECHELLE	TITRE / *TITLE*		
JMA	10/02/2022	1:1	Embase_PRM-FOEM		
MATERIAU / *MATERIAL*		SYSTEME	REFERENCE	FINITION / *FINISH*	COULEUR / *COLOR*
			M_FLG_01269		


O'kabé bureaux, 67 avenue de Fontainebleau
94270 Le Kremlin-Bicêtre | France


FLUIGENT
SMART MICROFLUIDICS™

TOLERANCES / *TOLERANCES*
+/- 0.5 mm

FORMAT	LANGUE	PROPRIETE		PAGE
A4		All proprietary rights reserved by Fluigent. Unauthorized reproduction will be legally pursued.		1/3







-	-	Voir Feuille1	-
Indice	Date	Modification	Approuvé

Table de Révision

CREE / *CREATED*	DATE	ECHELLE	TITRE / *TITLE*		
JMA	10/02/2022	1:1	Embase_PRM-FOEM		
MATERIAU / *MATERIAL*		SYSTEME	REFERENCE	FINITION / *FINISH*	COULEUR / *COLOR*
			M_FLG_01269		

O'kabé bureaux, 67 avenue de Fontainebleau
94270 Le Kremlin-Bicêtre | France


FLUIGENT
SMART MICROFLUIDICS™

TOLERANCES / *TOLERANCES*
+/- 0.5 mm

FORMAT	LANGUE	PROPRIETE		PAGE
A4		All proprietary rights reserved by Fluigent. Unauthorized reproduction will be legally pursued.		2/3

PRESSURE CONNECTORS

- Push-Pull configuration:
Push-in connector
M5 thread - 4mm OD tubing
- FlowEZ configuration:
empty



Barb connector
M5 thread - 3/32" ID tubing

Push-in connector
M5 thread - 4mm OD tubing

Push-in connector
M5 thread - 4mm OD tubing

Indice	Date	Modification	Approuvé
-	-	Voir Feuille1	-

Table de Révision

CREE / *CREATED*	DATE	ECHELLE	TITRE / *TITLE*		
JMA	10/02/2022	1:1	Embase_PRM-FOEM		
MATERIAU / *MATERIAL*		SYSTEME	REFERENCE	FINITION / *FINISH*	COULEUR / *COLOR*
			M_FLG_01269		

O'kabé bureaux, 67 avenue de Fontainebleau
94270 Le Kremlin-Bicêtre | France

FLUIGENT
SMART MICROFLUIDICS

TOLERANCES / *TOLERANCES*				
+/- 0.5 mm				
FORMAT	LANGUE	PROPRIETE		PAGE
A4		All proprietary rights reserved by Fluigent. Unauthorized reproduction will be legally pursued.		3/3